BioStrip HCG

One Step Urine Pregnancy test

Catalogue# 4220-5

INTENDED USE

For the rapid determination of human chorionic gonadotropin (hCG) in urine specimens. This test kit is used to obtain a visual, qualitative result for early detection of pregnancy. This product is for professional use only.

SUMMARY AND EXPLANATION

Human chorionic gonadotropin (hCG) is a glycoprotein hormone secreted by the developing placenta shortly after implantation. hCG can be detected in the urine and serum of pregnant women as early as 6 to 15 days after conception. The concentration of hCG increases to 5-50 mIU/mL one week post implantation, reaches about 100 mIU/mL at the time of the first missed menstrual period, and peaks at 100,000-200,000 mIU/mL at the end of the first trimester¹.

The appearance of hCG soon after conception and its subsequent rise in concentration during early gestation growth make it an excellent marker for the early detection of pregnancy.

The BioStrip HCG test kit is a rapid test to detect the presence of hCG in urine specimens in a qualitative format sensitive to 20 mIU hCG/mL. The test utilizes a combination of monoclonal and polyclonal antibody reagents to selectively detect elevated level of hCG in urine. The immunological specificity of the test kit virtually eliminates cross reactivity interferences from the structurally related glycoprotein hormones hFSH, hLH, hTSH at physiological levels.

PRINCIPLE

The BioStrip HCG test is a qualitative, two site sandwich immunoassay for the determination of hCG in urine. The membrane was precoated with goat anti-hCG on the test band region and goat anti-mouse on the control band region. During the test, the patient urine is allowed to react with the coloured colloidal gold-monoclonal anti-hCG conjugate which was predried on the test strip. The mixture then moves upward on the membrane by the capillary action. For a positive result, a colour band with the specific antibody-hCG-coloured conjugate complex will form on the membrane. Absence of this coloured band in the test band region suggest a negative result. To serve as a procedural control, a coloured band at control region will always appear regardless the presence of hCG in test specimen.

STORAGE AND STABILITY

The test kit can be stored at 2°C to 30°C in the desiccated container up to the expiration date. The test strips should be kept from direct sunlight, moisture and heat.

The cap of the container houses a desiccant to keep the test strips dry. Replace cap tightly, immediately after use.

CAUTION

1. For professional use only
2. For in vitro diagnostic use only.
3. Do not use test kit beyond expiry date.
4. For Single use only, do not reuse

REAGENTS AND MATERIALS SUPPLIED

Test strips: Strip coated with goat anti-hCG and coloured colloidal gold - monoclonal mouse anti-hCG conjugated predried in pad.

SPECIMEN COLLECTION


Any urine specimen is appropriate for hCG testing, but the first morning urine is optimal because it generally contains the highest concentration of hCG. Urine specimens may be collected in any clean, dry, plastic or glass container. If specimens cannot be assayed immediately, they may be stored in the refrigerator at 2°C to 8°C for up to 72 hours prior to assay. Preservatives are not required.

Specimen should be allowed to reach room temperature before testing. Very turbid urine specimen should be centrifuged or filtered prior to use. Samples with visible precipitates should be allowed to settle and the clear supernatant used for testing.

TEST PROCEDURE


Test strips, patient samples, controls and reference materials should be brought to room temperature 20°C to 30°C prior to testing.

1. Remove the test strip from the desiccated container. Cap the container tightly and promptly immediately after use. Label the strip with patient or control identification.
2. Immerse the strip into the urine with the arrow end pointing towards the urine. Do not immerse pass the MAX (maximum) line. You may leave the strip in the urine or you may take the strip out after a minimum of 10 seconds and lay the strip flat on a clean, dry, non-absorbent surface (e.g. mouth of the urine container).
3. Coloured band(s) will appear. Depending on the concentration of hCG in the test specimen, positive results may be observed in as short as 10-30 seconds. However, confirm negative results, the complete reaction time of 5 minutes is required. **Do not interpret results after 10 minutes.**

INTERPRETATION OF RESULTS

1. **Positive:** Distinct colour bands appear on the control and test region. Colour intensity of the test bands may vary.
2. **Negative:** Only one colour band appears on the control region. No apparent band on the test region.
3. **Invalid:** No visible band at all or no visible band at the control region. Repeat test.

QUALITY CONTROL

The procedural control is included in the test. A colour band appearing on the control region of the strip indicates proper performance and reactive reagents.

Good laboratory practices include the use of control specimens to ensure proper kit performance. Negative and positive controls containing hCG at various concentrations are available commercially.

PROCEDURE LIMITATION

1. A number of conditions other than pregnancy including trophoblastic disease and certain nontrophoblastic neoplasms cause elevated levels of hCG. These diagnosis should be considered if appropriate to the clinical evidence.
2. If the urine specimen is too dilute (i.e. low specific gravity) it may not contain representative level of hCG. If pregnancy is still suspected, a first morning urine should be obtained from the patient 48-72 hours later and tested.
3. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after clinical and laboratory findings have been evaluated.
4. False negative reaction may occur when the HCG level in the urine specimen is extremely high. This typically occurs when the patient's HCG levels peaks at 13-16 week pregnancy. If this is suspected, one can dilute one drop of urine with 10 drops of saline and retest the diluted urine using a fresh test strip. If the diluted specimen gives a positive result, one can conclude the pregnancy test is positive.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity:

The BioStrip HCG detects urine hCG concentrations greater than

20 mIU/mL as indicated by the appearance of a colour band at the test region.

Specificity of the BioStrip HCG test was determined from cross reaction studies with known amounts of luteinizing hormone (hLH), follicle stimulating hormone (hFSH), and thyroid stimulating hormone (hTSH). 200 mIU hLH/mL, 1,000 mIU hFSH/mL and 1000 μ IU hTSH/mL all gave negative results.

Standardization

The BioStrip HCG test has been standardised to World Health Organization Third International Standard (WHO 3rd IS 75/589).

Interference Testing

The following substances were added in hCG free and 25 mIU hCG/mL spiked urine samples. None of the substances at concentration tested interfered in the assay.

Acetaminophen	20 mg/mL
Acetylsalicylic acid	20 mg/mL
Ascorbic acid	20 mg/mL
Atropine	20 mg/mL
Caffeine	20 mg/mL
Genesic acid	20 mg/mL
Glucose	2 g/dL
Haemoglobin	1 mg/dL

DISPOSAL

Urine samples and used test strips are potentially infectious. Follow each laboratory's established procedure for proper handling and disposal.

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Innovatek Medical Inc.