# QuickStep® Plus HCG Combo

## Serum/Urine Pregnancy test

Catalogue # 4220-100-27

### INTENDED USE

For the rapid determination of human chorionic gonadotropin (hCG) in urine or serum specimens. This test kit is used to obtain a visual, qualitative result for early detection of pregnancy. For professional use only.

### SUMMARY AND EXPLANATION

Human chorionic gonadotropin (hCG) is a glycoprotein hormone secreted by the developing placenta shortly after implantation. hCG can be detected in the urine and serum of pregnant women as early as 6 to 15 days after conception. The concentration of hCG increases to 5-50 IU/L one week post implantation, reaches about 100 IU/L at the time of the first missed menstrual period, and peaks at 100,000-200,000 IU/L at the end of the first trimester<sup>1</sup>.

The appearance of hCG soon after conception and its subsequent rise in concentration during early gestation growth make it an excellent marker for the early detection of pregnancy.

The QuickStep® Plus HCG Combo test kit is a rapid test to detect the presence of hCG in urine or serum specimens in a qualitative format sensitive to 20 mIU hCG/mL. The test utilizes a combination of monoclonal and polyclonal antibody reagents to selectively detect elevated level of hCG in urine or serum. The immunological specificity of the test kit virtually eliminates cross reactivity interferences from the structurally related glycoprotein hormones hFSH, hLH, hTSH at physiological levels.

### TEST PRINCIPLE

The QuickStep® Plus HCG Combo test is a qualitative; two site sandwich immunoassay for the determination of hCG in urine or serum. The membrane was pre-coated with goat anti-hCG on the test band region and goat anti-mouse on the control band region. During the test, the patient urine is allowed to react with the coloured colloidal gold-mono-clonal anti-hCG conjugate which was pre-dried on the test strip. The mixture then moves upward on the membrane by the capillary action. For a positive result, a colour band with the specific antibody-hCG-coloured conjugate complex will form on the membrane. Absence of this coloured band in the test band region suggests a negative result. To serve as a procedural control, a coloured band at control region will always appear regardless the presence of hCG in test specimen.

### STORAGE AND STABILITY

The test kit can be stored at room temperature from +2° to 30°C in the sealed pouch up to the expiration date. The test should be kept from direct sunlight, moisture and heat.

### CAUTION

1. For professional use only.
2. For in vitro diagnostic use only.
3. Do not use test kit beyond expiry date.
4. For single use only, do not reuse.

### REAGENTS AND MATERIALS SUPPLIED

QuickStep Plus HCG Combo Serum/Urine Pregnancy test device: Containing membrane coated with goat anti-hCG and coloured

colloidal gold - monoclonal mouse anti-hCG conjugated predried in pad.


### SPECIMEN COLLECTION AND HANDLING

Any urine or serum specimen is appropriate for hCG testing, but the first morning urine is optimal because it generally contains the highest concentration of hCG. Urine specimens may be collected in any clean, dry, plastic or glass container. Collect blood and allow to clot, separate serum for testing. If specimens cannot be assayed immediately, they may be stored in the refrigerator at 2° to 8°C for up to 72 hours prior to assay. Preservatives are not required.

Specimen should be allowed to reach room temperature before testing. Very turbid urine specimen should be centrifuged or filtered prior to use. Samples with visible precipitates should be allowed to settle and the clear supernatant used for testing. Very lipemic serum and plasma is not suitable for testing.

### TEST PROCEDURE

1. Review Specimen Collection instructions. Test device, patient's samples, or controls should be brought to room temperature (18° to 25°C) prior to testing.
2. Remove the test device from its protective pouch. Label the device with patient or control identifications.
3. Add 3 drops of urine or serum into the samples well. Hold the pipette in a straight up and down position - not at an angle. For each sample or control, use a separate pipette and test device.


Add specimen

4. Wait for coloured bands to appear. Depending on the concentration of hCG, positive results may be observed in as short as 15 seconds. However, to confirm negative results, the complete reaction time (5 minutes for urine, 10 minutes for serum) is required. **Do not interpret results after 20 minutes.**


### INTERPRETATION OF RESULTS

**Negative:** Only one colour band appears on the Control (C) region. No apparent band on the Test (T) region.


Negative

**Positive:** Distinct colour bands appear on the Control (C) and Test (T) region. Colour intensity of the bands may vary.

Positive

**Invalid:** No visible band at all or no visible band at the Control

(C) region. Repeat test with a fresh test device.

Invalid

### QUALITY CONTROL

The procedural control is included in the test. A colour band appearing on the control region of the strip indicates proper performance and reactive reagents.

Good laboratory practices include the use of control specimens to ensure proper kit performance. Negative and positive controls containing hCG at various concentrations are available commercially.

### LIMITATIONS

1. A number of conditions other than pregnancy including trophoblastic disease and certain non-trophoblastic neoplasms cause elevated levels of hCG. These diagnosis should be considered if appropriate to the clinical evidence.
2. Patient with high level natural antibodies against mouse and those receiving therapeutic monoclonal antibodies of mouse origin may give false positive test results.
3. If the urine specimen is too dilute (i.e. low specific gravity) it may not contain representative level of hCG. If pregnancy is still suspected, a first morning urine should be obtained from the patient 48-72 hours later and tested.
4. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after clinical and laboratory findings have been evaluated.
5. At very high level of hCG (end of first trimester), Hook Effect may cause the test band to become very weak or false negative result. The urine specimen can be diluted 1:100 with saline and tested again using a fresh test device. The diluted specimen should give a strong positive test result.

### PERFORMACNE CHARACTERISTICS

#### Sensitivity and Specificity:

The QuickStep® Plus HCG Combo detects urine or serum hCG concentrations greater than 20 IU/L as indicated by the appearance of a colour band at the test region.

Specificity of the QuickStep® Plus test was determined from cross reaction studies with known amounts of luteinizing hormone (hLH), follicle stimulating hormone (hFSH), and thyroid

stimulating hormone (hTSH). 300 IU hLH/L, 1,000 IU hFSH/L and 1,000 µIU hTSH/mL all gave negative results.

When tested with 200 known negative sera, one sample gave very weak positive result.

### Standardisation

The Quickstep® Plus test has been standardised to World Health Organisation First International Reference Preparation (WHO 1st IRP 75/537).

### Interference Testing

The following substances were added in hCG free and 50IU hCG/L spiked urine samples. None of the substances at concentration tested interfered in the assay.

Acetaminophen	20 mg/mL
Acetylsalicylic acid	20 mg/mL
Ascorbic acid	20 mg/mL
Atropine	20 mg/mL
Caffeine	20 mg/mL
Gentesic acid	20 mg/mL
Glucose	2 g/dL
Haemoglobin	1 mg/dL

### DISPOSAL

Urine, serum samples and used test devices are potentially infectious. Follow each laboratory's established procedure for proper handling and disposal.

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**Innovatek Medical Inc.**