QuickStick HCG Combo

Serum / Urine Pregnancy Test

Catalogue #4220-100-24

INTENDED USE

For the rapid qualitative determination of human chorionic gonadotropin (hCG) in urine and serum specimens. This test kit is used to obtain a visual, qualitative result for early detection of pregnancy. This product is for professional use only.

SUMMARY AND EXPLANATION

Human chorionic gonadotropin (hCG) is a glycoprotein hormone secreted by the developing placenta shortly after implantation. hCG can be detected in the urine and serum of pregnant women as early as 6 to 15 days after conception. The concentration of hCG increases to 5-50 mIU/mL one week post implantation, reaches about 100 mIU/mL at the time of the first missed menstrual period, and peaks at 100,000-200,000 mIU/mL at the end of the first trimester.

The appearance of hCG soon after conception and its subsequent rise in concentration during early gestation growth make it an excellent marker for the early detection of pregnancy.

The QuickStick HCG Combo test kit is a rapid test to detect the presence of hCG in serum/urine specimens in a qualitative format sensitive to 20 mIU hCG/mL. The test utilizes a combination of monoclonal and polyclonal antibody reagents to selectively detect elevated level of hCG in serum/urine. The immunological specificity of the test kit virtually eliminates cross reactivity interferences from the structurally related glycoprotein hormones hFSH, hLH, hTSH at physiological levels.

PRINCIPLE

The QuickStick HCG Combo test is a qualitative, two site sandwich immunoassay for the determination of hCG in serum/urine. The membrane was precoated with goat anti-hCG on the test band region and goat anti-mouse on the control band region. During the test, the patient serum/urine is allowed to react with the coloured colloidal gold-monoclonal anti-hCG conjugate which was predried on the test strip. The mixture then moves upward on the membrane by the capillary action. For a positive result, a colour band with the specific antibody-hCG-coloured conjugate complex will form on the membrane. Absence of this coloured band in the test band region suggest a negative result. To serve as a procedural control, a coloured band at control region will always appear regardless the presence of hCG in test specimen.

STORAGE AND STABILITY

The test kit can be stored at 2° to 3°C in the sealed pouch up to the expiration date. The test strips should be kept from direct sunlight, moisture and heat.

WARNING AND PRECAUTION

1. For Professional use only
2. For in vitro diagnostic use only.

3. Do not use test kit beyond expiry date.
4. For single use only, do not reuse

REAGENTS AND MATERIALS SUPPLIED

Test strips: Strip coated with goat anti-hCG and coloured colloidal gold - monoclonal mouse anti-hCG conjugated predried in pad

SPECIMEN COLLECTION

Any serum/urine specimen is appropriate for hCG testing, but the first morning urine is optimal because it generally contains the highest concentration of hCG. Urine specimens may be collected in any clean, dry, plastic or glass container.

For serum, collect blood without anticoagulant. Allow the blood to clot and collect the serum. Lipemic and severely haemolysed serum is not suitable. Do not use plasma for testing.


If specimens cannot be assayed immediately, they may be stored in the refrigerator at 2° to 8°C for up to 72 hours prior to assay. Preservatives are not required.

Specimen should be allowed to reach room temperature before testing. Very turbid urine specimen should be centrifuged or filtered prior to use. Samples with visible precipitates should be allowed to settle and the clear supernatant used for testing.

TEST PROCEDURE


Test strips, patient samples, controls and reference materials should be brought to room temperature 20°C to 30°C prior to testing.

1. Remove the test strip from the sealed pouch. Label the strip with patient or control identification.
2. Immerse the strip into the urine or serum specimen with the arrow end pointing towards the urine. Do not immerse pass the MAX (maximum) line. You may leave the strip in the specimen or you may take the strip out after a minimum of 10 seconds and lay the strip flat on a clean, dry, non-absorbent surface.
3. Coloured band(s) will appear. Depending on the concentration of hCG in the test specimen, positive results may be observed in as short as 10-30 seconds. However, to confirm negative results, the complete reaction time of 5 minutes is required. Do not interpret results after 10 minutes.

INTERPRETATION OF RESULTS

1. **Positive:** Distinct colour bands appear on the control and test region. Colour intensity of the test bands may vary.
2. **Negative:** Only one colour band appears on the control region. No apparent band on the test region.
3. **Invalid:** No visible band at all or no visible band at the control region. Repeat test.

QUALITY CONTROL

The procedural control is included in the test. A colour band appearing on the control region of the strip indicates proper performance and reactive reagents.

Good laboratory practices include the use of control specimens to ensure proper kit performance. Negative and positive controls containing hCG at various concentrations are available commercially.

PROCEDURE LIMITATION

1. A number of conditions other than pregnancy including trophoblastic disease and certain nontrophoblastic neoplasms cause elevated levels of hCG. These diagnosis should be considered if appropriate to the clinical evidence.
2. If the urine specimen is too dilute (i.e. low specific gravity) it may not contain representative level of hCG. If pregnancy is still suspected, a first morning urine should be obtained from the patient 48-72 hours later and tested.
3. As with all diagnostic tests, a definitive clinical diagnosis should not be based on the results of a single test, but should only be made by the physician after clinical and laboratory findings have been evaluated.
4. False negative reaction may occur when the HCG level in the urine specimen is extremely high. This typically occurs when the patient's HCG levels peaks at 13-16 week pregnancy. If this is suspected, one can dilute one drop of urine with 10 drops of saline and retest the diluted urine using a fresh test strip. If the diluted specimen gives a positive result, one can conclude the pregnancy test is positive.

PERFORMANCE CHARACTERISTICS

Sensitivity and Specificity:

The QuickStick HCG Combo detects urine/serum hCG concentrations greater than 20 mIU/mL as indicated by the appearance of a colour band at the test region.

Specificity of the QuickStick HCG Combo test was determined from cross reaction studies with known amounts of luteinizing hormone (hLH), follicle stimulating hormone (hFSH), and thyroid stimulating hormone (hTSH). 200 mIU hLH/mL, 1,000 mIU hFSH/mL and 1,000 mIU hTSH/mL all gave negative results. When tested with 200 known negative sera samples, one sample gave very weak positive result.

Standardization

The QuickStick Combo test has been standardised to World Health Organization Third International Standard (WHO 3rd IS 75/589).

Interference Testing

The following substances were added in hCG free and 25 mIU hCG/mL spiked urine samples. None of the substances at concentration tested interfered in the assay.

Acetaminophen	20 mg/mL
Acetylsalicylic acid	20 mg/mL
Ascorbic acid	20 mg/mL
Atropine	20 mg/mL
Caffeine	20 mg/mL
Genesic acid	20 mg/mL
Glucose	2 g/dL
Haemoglobin	1 mg/dL

DISPOSAL

Urine / Serum samples and used test strips are potentially infectious. Follow each laboratory's established procedure for proper handling and disposal.

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Innovatek Medical Inc.