AMNICATOR® Test

For detecting ruptured amniotic membrane

INTENDED USE

For in vitro diagnostic use. Amnicator® is a rapid and simple method intended for detecting ruptured amniotic membrane in pregnant women. To be used by qualified healthcare professionals.

SUMMARY AND EXPLANATION

Rupture of the amniotic membrane can result in the leaking of amniotic fluid into the vagina. The presence of amniotic fluid tends to elevate the pH of the upper vagina. Detection of this pH increase using a pH indicator dye has been shown to assist in determining the presence of amniotic fluid as an indirect indication of ruptured amniotic membrane.

PRINCIPLE OF THE PROCEDURE

Amnicator® is individually wrapped sterile swabs impregnated with a sensitive pH indicator, nitrazine yellow. The swab absorbs fluid associated with the tissue and the dye develops a colour which correlates with the pH of the absorbed fluid.

The principle of the Amnicator® test is that normal vaginal pH at term ranges from 4.5 to 6.0 and amniotic fluid is usually 7.0 to 7.5. By means of a simple pH colour change of nitrazine yellow, which has a sharp end point (at pH 6.4 to 6.8) exactly in the desired pH range, one can accurately detect the presence of amniotic fluid in the vagina.

STORAGE AND STABILITY

Amnicator® test in individually wrapped pouches should be stored at room temperature (15° to 30°C) in the box. Do not expose the swabs to bright light. Product stored under these conditions will be stable until the expiry date shown on the label. Do not open the pouch until ready to use.

PRECAUTIONS

1. Amnicator® is intended for in vitro use only.
2. Do not open the pouch until ready to use.
3. Do not use Amnicator® after expiry date shown on the product label.
4. Following contact with the vagina, swabs should be considered potentially infectious and precautions appropriate for bio-hazards must be observed.
5. For single use only. Do not reuse swabs.
6. The swab is sterile in the unopened and undamaged pouch.
7. The procedures, storage conditions, precautions and limitations specified in these directions must be adhered to in order to obtain valid information.

REAGENT AND MATERIALS SUPPLIED


Amnicator® sterile swab in an individually wrapped peel pouch.

TEST PROCEDURE

1. Cleanse the vulva and vaginal introitus.
2. Carefully insert a bi-value speculum (e.g. Cusco's speculum) into the vagina exposing the cervix.
3. Pass an Amnicator® into the vagina and sample for 10 to 15 seconds any liquid present at the external cervical os or in posterior fornix of the vagina.
4. Withdraw the Amnicator® swab and inspect the colour of the swab tip immediately.
5. Compare the colour of the Amnicator® swab with the colour chart.

INTERPRETATION OF RESULTS

Compare the colour of the Amnicator® swab tip with this colour chart

QUALITY CONTROL

Laboratories wishing to perform optional in-house quality control may employ buffers corresponding to the pH values listed on the Amnicator® colour chart. Remove one Amnicator® swab from its sleeve. Wet the swab's tip with 3 to 4 drops of buffer solution. Immediately compare the colour developed on the swab tip to the closest matching colour on the Amnicator® colour chart. If the pH value written next to the colour selected on the card corresponds to the pH of the buffer solution then the swab is performing as expected. If the pH indicated on the Colour card fails to match the pH of the buffer used, repeat the test with a fresh swab. If the pHs still do not match, the swab is not performing as expected and the remaining kit should not be used to test clinical specimens. Individual swabs tested in this manner should be discarded and must not be used for testing clinical samples.

LIMITATIONS

1. Amnicator® test is designed to be used by qualified medical professional and is intended as an aid to professional diagnosis.
2. Amnicator® test can only indicate a change in pH value and should be used only as indicated in the Test Procedure described above.
3. Antibiotic therapy or presence of vaginal infection can alter the normal vaginal pH which may lead to false positive result.
4. If vaginal infection is suspected, microbiological testing should be used to exclude infection.

REFERENCES

1. Abe T. The detection of the rupture of foetal membranes with the nitrazine indicator. *Am J. Obst. & Gynec.* 1940. 39:400-404.
2. Drife J. Preterm rupture of the membranes. *Brit. Med. J.* 1982; 285:583.
3. Mills A., Ganloch D. Use of the nitrazine yellow swab test in the diagnosis of the ruptured membranes. *Brit. J. Obst. & Gynec.* 1977; 84:138-140.
4. Pauersterin C. Premature rupture of the membranes. In: *Clinical Obstetrics*, Anonymous, ed. John Wiley & Sons and Churchill Livingstone. 1987; 367-381.
5. Pritchard J, Macdonald P. *Williams Obstetrics*, anonymous, ed. New York: Appleton-Century-Crofts, 1960; 407-408.
6. Filet JP, More N, Librati C, Ruffie A, Delouis P, Cluzeau MH, Hocke C, and Lang JJ. Évaluation de trois méthodes diagnostiques dans la rupture prématrice des membranes. *Rev. Fr. Gynecol. Obstet.* 1994; 89:123-128.
7. Dickinson, F., Soltani H. & Allsop J., 2005, Review of the Use of Nitrazine Yellow Swabs in Pre-Labour. 27th Triennial Congress International Confederation of Midwives, Brisbane, Australia.

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