



Pre-moistened sponge swab for microbial surface sampling, an efficient, ready-to-use surface sampling solution.

The POLYWIPE™ sponge swab is a sterile, pre-moistened blue sponge designed for the microbial surface sampling of a wide variety of surfaces in industries such as food manufacturing, and infection prevention & control.



Rapid sampling allows work to continue at any location and transport to the lab.



Reducing cost per test while managing waste for a positive sustainable impact.



Tamper proof packaging minimises expensive retesting. Different specimens require different solutions for optimum results.

## TRADITIONAL FORMATS

### PEEL POUCH & CONTAINER

-  Ready-to-use and individually sealed.
-  Suitable for quick, efficient sampling.
-  Multiple buffer options available.



## NEW FORMAT

### EASY OPEN STOMACHER BAG

The Easy Open Stomacher Bag and POLYWIPE™ sponge swab offers sample security while providing enhanced convenience through reduced plastics and increased sampling efficiency.



Easy Open Stomacher Bag doubles as a glove for easier handling during sampling with less packaging waste for a positive environmental impact.



The same Stomacher Bag is used for sample processing, improving efficiency and waste reduction.



Lightweight pre-moistened samples are easily transported from samplings locations to the labs.



The Stomacher Bag has a large area designated for writing sample details. Removing the need for additional labels or paperwork.



# POLYWIFE™

ENVIRONMENTAL SAMPLING



## WHAT YOU NEED TO KNOW



- Tamper Proof Seal
- Shelf-Life: 12 months
- Stomacher Bag & Template Pack Size: 40
- Peel Pouch Pack Size: 50
- Container Pack Size: 30

### Sponge Dimensions:

W: 45mm  
H: 98mm  
D: 4mm

### Format A: Grip Seal Bags

Format B:  
Grip Seal Bags &  
Disposable Gloves

### Format TB: Twist Bags

| Most Common Use                       | Product Code | Blue Sponge Swab Description  |
|---------------------------------------|--------------|---|
| Multipurpose                          | M451001      | Phosphate Buffer - Easy Open Stomacher Bag  |
| Neutralising Cleaning Residue         | M451002      | NRS+ - Easy Open Stomacher Bag  |
| Saline Non-Neutralising               | M451003      | 0.9% Saline - Easy Open Stomacher Bag   |
| Carcass                               | M451004      | Peptone Saline MRD - Easy Open Stomacher Bag  |
| Salmonella Testing                    | M451005      | Buffered Peptone Water - Easy Open Stomacher Bag  |
| Flat Surfaces                         | M451101      | Sampling Template - 10x10cm / 100 sqcm - Individually Wrapped                                     |
| Salmonella Testing                    | MW722        | Buffered Peptone Water - Peel Pouch   |
| Salmonella Testing                    | MW722A       | Buffered Peptone Water - 50 Grip Seal Bags 210x130mm - Peel Pouch                                 |
| Salmonella Testing                    | MW722B       | Buffered Peptone Water - 50 Grip Seal Bags 210x130mm - 50 Disposable Gloves - Peel Pouch          |
| Neutralising Residue                  | MW723A       | Neutralising Rinse Solution NRS - 50 Grip Seal Bags 210x130mm - Peel Pouch                        |
| Neutralising Residue                  | MW723B       | Neutralising Rinse Solution NRS - 50 Grip Seal Bags 210x130mm - 50 Disposable Gloves - Peel Pouch |
| Neutralisation of Quaternary Ammonium | MW724        | Lethen Broth - Peel Pouch   |
| Neutralisation of Quaternary Ammonium | MW724A       | Lethen Broth - 50 Grip Seal Bags 210x130mm - Peel Pouch   |
| Neutralisation of Quaternary Ammonium | MW724TB      | Lethen Broth - 50 Twist Bags - Peel Pouch   |
| Multipurpose                          | MW725        | Phosphate Buffer - 110ml Polycon Container Blue Cap   |
| Carcass                               | MW726        | Peptone Saline MRD - Peel Pouch   |
| Carcass                               | MW726A       | Peptone Saline MRD - 50 Grip Seal Bags 210x130mm - Peel Pouch                                     |
| Carcass                               | MW726B       | Peptone Saline MRD - 50 Grip Seal Bags - 50 Disposable Gloves - Peel Pouch                        |
| Multipurpose                          | MW728        | Phosphate Buffer - 50 Disposable Gloves - 110ml Polycon Container Blue Cap                        |
| Multipurpose                          | MW729        | Phosphate Buffer - Peel Pouch   |
| Multipurpose                          | MW729A       | Phosphate Buffer - 50 Grip Seal Bags 210x130mm - Peel Pouch                                       |
| Multipurpose                          | MW729B       | Phosphate Buffer - 50 Grip Seal Bags 210x130mm - 50 Disposable Gloves - Peel Pouch                |
| Multipurpose                          | MW729C       | Phosphate Buffer - 100 Disposable Gloves - Peel Pouch   |
| Multipurpose                          | MW729TB      | Blue Sponge Swab - Phosphate Buffer - 50 Twist Bags - 100 Disposable Gloves - Peel Pouch          |



We ensure every sample is collected and transported in optimal condition, delivering results that advance care, discovery, and progress.

## HAVE YOU THOUGHT ABOUT:

The STERISWAB™ for sampling critical surfaces including nozzles & equipment in sterile units in hospitals, manufacturing, and pharmaceutical industries.